

Serological tests 2
(Antigen antibody interactions)
Lab 3


Samer Alqaraleh,
PhD. Nanobiotechnology
Faculty of Medicine, Mutah university
Immunology, 2nd year students

HCG Latex Agglutination Test



Lab. 6





DIAGNOSTIC AUTOMATION, INC.
IMMUNODIAGNOSTICS



MICROWELL ELISA 
Brucella IgG

REF 1501-8

96 Test

LOT BRG-196

80°

2011-12

LOT

DIAGNOSTIC AUTOMATION, INC.
10000 WILSON BLVD, SUITE 100, FORT WORTH, TX 76154

USA

10



CODE
J79

PACK
25T

RHEUMATOID FACTOR

LATEX TEST KIT (R.A. TEST KIT)

Contents

Reagent 1 : RF Antigen (Gamma Globulin)

Reagent 2 : Positive Control

Reagent 3 : Negative Control

Accessories : Disposable Plastic Dropper, Disposable Applicator
Sticks, Rubber Teat, Glass Slides

Shake well before use



ISO 9001 & 13485
Certified Company



Mfg. Lic. No. G-639



BEACON
DIAGNOSTICS PVT. LTD.

424, NEW GIDC, KASBPUR, NAVSARI - 395 404, GUJ. INDIA
Tel : +91 - 2637 - 202040 | Fax : www.beaconindia.com

ANTI - A, ANTI - B, ANTI - A


METHOD - IMMUNOAGGLUTINATION
PRODUCT CODE - AB10 | AB11

INSTRUCT

INTENDED USE: Test for qualitative detection of ABO ar


INTRODUCTION

The mouse monoclonal IgM Anti-A, Anti-B and Anti-AB antibodies



ANAMOL
THE ORIGINAL MAKERS

Making Biochemical Formulations Since 1981



Precipitation Reaction Types

They are of mainly three types:

1. Precipitation in Solution

- a) Ring Test.
- b) Slide Test.
- c) Tube Test


2. Precipitation in Agar.

- a) Single diffusion immunodiffusion test (Mancini test)
- b) Double diffusion immunodiffusion test (Ouchterlony test)


3. Precipitation in Agar in an electric field (immunoelectrophoresis).

Single


Single radial immunodiffusion



1



$$r \propto [Ag]$$




Method

- Ab in gel
- Ag in a well
- Quantitative
- Interpretation



1. Diameter of ring is proportional to the concentration

2. Height of the Rocket




Electroimmunodiffusion

Double



Ouchterlony Gel Diffusion

- Both antigen and antibody can diffuse independently.
- Holes punched in agar.
- Known antibody or antigen added to center well.
- Known sample added to outer well.
- Unknown sample added to outer well next to unknown sample.
- Wait for bands to form.




Principle of Immunoelectrophoresis

- When an electric current is applied to a slide layered with gel, the antigen mixture placed in wells is separated into individual antigen components according to their charge and size.
- Following electrophoresis, the separated antigens are reacted with specific antisera placed in troughs parallel to the electrophoretic migration, and diffusion is allowed to occur. Antiserum present in the trough moves toward the antigen components resulting in the formation of separate precipitin lines in 18-24 hrs, each indicating reaction between individual proteins with their antibody.

Immuno-electrophoresis

Gel Tank (with Agarose Gel)

Power Supply






FIG. 14-6. Counter-current immunoelectrophoresis.




FIG. 14-7. Rocket electrophoresis.

Neutralization Tests

Neutralization is an antigen–antibody reaction in which the biological effects of viruses and toxins are neutralized by homologous antibodies known as neutralizing antibodies. These tests are broadly of two types:

- (a) virus neutralization tests
- (b) toxin neutralization tests.

Virus neutralization tests

- Neutralization of viruses by their specific antibodies are called virus neutralization tests. Inoculation of viruses in cell cultures, eggs, and animals results in the replication and growth of viruses. When virus-specific neutralizing antibodies are injected into these systems, replication and growth of viruses is inhibited.
- **Viral hemagglutination inhibition** test is an example of virus neutralization test frequently used in the diagnosis of viral infections, such as influenza, mumps, and measles. If patient's serum contains antibodies against certain viruses that have the property of agglutinating the red blood cells, these antibodies react with the viruses and inhibit the agglutination of the red blood cells.

Toxin neutralization tests

- Toxin neutralization tests are based on the principle that biological action of toxin is neutralized on reacting with specific neutralizing antibodies called antitoxins. Examples of neutralization tests include:
- *In vivo*— Schick test to demonstrate immunity against diphtheria and *Clostridium welchii* toxin neutralization test in guinea pig or mice.
- *In vitro*— (a) anti-streptolysin O test and (b) Nagler reaction used for rapid detection of *C. welchii*.

Immunofluorescence

The property of certain dyes absorbing light rays at one particular wavelength (ultraviolet light) and emitting them at a different wavelength (visible light) is known as fluorescence. Fluorescent dyes, such as fluorescein isothiocyanate and lissamine rhodamine, can be tagged with antibody molecules. They emit blue-green and orange-red fluorescence, respectively under ultraviolet (UV) rays in the fluorescence microscope. This forms the basis of the immunological test.

Immunofluorescence tests have wide applications in research and diagnostics.

These tests are broadly of two types:

1. Direct immunofluorescence test
2. Indirect immunofluorescence test

Direct immunofluorescence test is widely used for detection of bacteria, parasites, viruses, fungi, or other antigens in CSF, blood, stool, urine, tissues, and other specimens.




FIG. 14-13. Direct fluorescent antibody test.

Indirect immunofluorescence test

Indirect immunofluorescence is a two-stage process.

- **First stage**, a known antigen is fixed on a slide. Then the patient's serum to be tested is applied to the slide, followed by careful washing. If the patient's serum contains antibody against the antigen, it will combine with antigen on the slide.
- **Second stage**, the combination of antibody with antigen can be detected by addition of a fluorescent dye-labeled antibody to human IgG, which is examined by a fluorescence microscope.

Lab 4

Enzyme Immunoassays

- Enzyme immunoassays (EIAs) can be used for detection of either antigens or antibodies in serum and other body fluids of the patient.
- In EIA techniques, antigen or antibody labeled with enzymes are used. Alkaline phosphatase, horseradish peroxidase, and galactosidase are the enzymes used in the EIA tests.
- The commonly used EIAs are enzyme-linked immunosorbent assays (ELISAs).
- These assays involve the use of an immunosorbent specific to either the antigen or antibody. Following the antigen–antibody reaction, chromogenic substrate specific to the enzyme (peroxidase, alkaline phosphatase, etc.) is added.
- The reaction is detected by reading the optical density.
- Usually, a standard curve based on known concentrations of antigen or antibody is prepared from which the unknown quantities are calculated.

There are different types of ELISA

(a) indirect ELISA,

(b) sandwich ELISA,

(c) competitive ELISA

	Antigen
	Antibody
	Enzyme conjugated antibody
	Enzyme substrate
	Product




FIG. 14-14. Indirect ELISA test.

	Antigen
	Antibody
	Enzyme conjugated antibody
	Enzyme substrate
	Product




FIG. 14-15. Sandwich ELISA test.

Competitive ELISA

- Used for the estimation of antibodies present in a specimen, such as serum.
- Principle of the test is that two specific antibodies, one conjugated with enzyme and the other present in test serum (if serum is positive for antibodies), are used.
- Competition occurs between the two antibodies for the same antigen.
- Appearance of color indicates a negative test (absence of antibodies), while the absence of color indicates a positive test (presence of antibodies).

In this test, the microtiter wells are coated with HIV antigen.

- The sera to be tested are added to these wells and incubated at 37°C and then washed.
- If antibodies are present in the test serum, antigen–antibody reaction occurs.
- The antigen– antibody reaction is detected by adding enzyme-labeled-specific HIV antibodies.
- In a positive test, no antigen is left for these antibodies to act.
- Hence, the antibodies remain free and are washed away during the process of washing.
- When substrate is added, no enzyme is available to act on it.
- Therefore, positive result indicates no color reaction. In a negative test, in which no antibodies are present in the serum, antigen in the coated wells is available to combine with enzyme-conjugated antibodies and the enzyme acts on the substrate to produce color.