

## BRUCELLA Agglutination Tube Test

Case three: prozone

Tube10	Tube9	Tube8	Tube7	Tube6	Tube5	Tube4	Tube3	Tube2	Tube1
Control-	+Control	1/1280	1/640	1/320	1/160	1/80	1/40	1/20	1/10
-	+	-	-	+	+	-	-	-	-



## BRUCELLA Agglutination Tube Test

Case two: negative reaction

Tube10	Tube9	Tube8	Tube7	Tube6	Tube5	Tube4	Tube3	Tube2	Tube1
Control-	Control+	1/1280	1/640	1/320	1/160	1/80	1/40	1/20	1/10
-	+	-	-	-	-	-	-	-	-

# BRUCELLA Agglutination Tube Test

Case one: positive reaction

Tube10	Tube9	Tube8	Tube7	Tube6	Tube5	Tube4	Tube3	Tube2	Tube1
Control-	Control+	1/1280	1/640	1/320	1/160	1/80	1/40	1/20	1/10
-	+	-	-	-	-	+	+	+	+



# COMPLEMENT FIXATION TEST

## positive case INFECTION

There is No haemolysis in the first tubes, and haemolysis in the rest of the test tube + haemolysis in the Antigen, Serum and Complement Control tubes.

The titre is the last tube showing No haemolysis.

Tube 1	Tube 2	Tube 3	Tube 4	Tube 5	Tube RBCs control	Tube Antigen Control	Tube Serum control	Tube Complement Control
NH	NH	NH	NH	H	NH	H	H	H

H: Haemolysis.

NH: No Haemolysis.




## LATEX AGGLUTINATION


TEST	SPECIMEN	TYPE OF REACTION SENSITIVITY	ANTIGEN	ANTIBODY	PRINCIPLE	CAUSATIVE AGENT	DISEASE CASE
<b>HCG</b> qualitative	urine / serum	passive agglutination	HCG Ag	anti-HCG Ab on latex	HCG + anti-HCG on latex = agglu. Within 2 min.		pregnancy
<b>CRP</b> qualitative & semiquantitative	serum	passive agglutination	CRP-C reactive protein	anti-CRP on latex	CRP+anti-CRP on latex = agglu. Within 2 min.	many bacteria & viruses	acute stage of inflammation diseases
<b>RF</b> qualitative & semiquantitative	serum	passive agglutination	human IgG on latex	IgM or IgG in serum (RF)	RF+human IgG on latex = agglu. Within 2 min.	UNKNOWN	Rheumatoid arthritis
<b>ASO</b> qualitative & semiquantitative	serum	passive agglutination	Streptolysin O on latex	anti-Streptolysin O in serum	Streptolysin O on latex + anti-Streptolysin O in serum = agglu. Within 2 min.	Streptococci pyogenes group A	streptococcal infections rheumatic fever



**Nucleolar:** the nucleoli stain as multiple solid bodies within the nucleus.



**Centromere:** large speckles of finite number, reactive antigen segregate with condensed chromosomes in cells undergoing mitosis.





6-examine for specific fluorescence under fluorescence microscope at a magnification of 200x or greater.


**Interpretation of the results**

The nuclear staining patterns observed with human epithelial cells (HEp)-2 provided include homogeneous, peripheral (rim), speckled, nucleolar and Centromere.

**Homogeneous:** the entire nucleus fluoresces evenly with a diffuse staining pattern.



**Peripheral (rim):** the nuclear membrane stains most intensely with a decreasing staining intensity of the nucleoplasm toward the center of the nucleus.



**Speckled:** Discrete coarse to fine round speckles fluoresce throughout the nucleus.

3

## COMPLEMENT FIXATION TEST

Anti-complementary reaction: NO haemolysis in all test tubes + NO Haemolysis in the serum, RBCs control tube.

The value of serum control tube is to test:

- The inability of serum alone to bind complement.
- The anti-complementary action of the serum, which is due to one of the followings:
  - a-presence of immune complexes in serum that binds to complement and prevent its haemolytic activity
  - b-Heparin therapy which inactivates complements due to consumption of Ca, Mg.
  - c-Old or contaminated serum has destructive action on complement.

Tube1	Tube2	Tube3	Tube4	Tube5	Tube RBCs control	Tube Ag control	Tube Serum control	Tube Complement control
NH	NH	NH	NH	NH	NH	H	NH	H



## COMPLEMENT FIXATION TEST

No infection (negative case): when patients serum does not contain antibodies, the complement will not be used and will be available to fix To and haemolyze the antibody coated red cells.

Haemolysis in all test tubes + haemolysis in the Ag , Ab,  
And complement control tube.

Tube1	Tube2	Tube3	Tube4	Tube5	Tube RBCs control	Tube Ag control	Tube Serum control	Tube Complement control
H	H	H	H	H	NH	H	H	H