

Microbiology

Lab 6

Antimicrobial susceptibility Test





General Microbiology Lab
Antimicrobial Susceptibility Test
Lab 6
2021-2022

Dr. Mohammad Odaibat
Department of Microbiology and Pathology
Faculty of Medicine, Mutah University



Medical Application



New antibiotics are continuously being developed



different bacteria acquire new resistant genes to the available antibiotics



determine the antibiotic susceptibility or resistance is required to determine most suitable antibiotic therapy

Methods of Antimicrobial Susceptibility Testing

1. Standardized filter-paper disc-agar diffusion (Kirby-Bauer method)



Qualitative Antimicrobial Susceptibility Testing

2. Minimum Inhibitory concentration (MIC)
& Minimum lethal concentration (MLC)

3. Epsilometer test (E-test)

**Quantitative Antimicrobial
Susceptibility Testing**

Standardized filter-paper disc-agar diffusion

Procedure




Urine sample



MacConkey agar



Gram negative bacilli
Lactose fermenter



Biochemical reactions

Antibiotic
susceptibility
test



E. coli




Glucose	A, G
Lactose	A, G
Maltose	A, G
Mannitol	A, G
Sucrose	A, G

indole	ve+
MR	ve+
VP	ve-
Citrate	ve-
Urease	ve-
H2S	ve-

Standardized filter-paper disc-agar diffusion

Principle



Mueller Hinton agar



Confluent growth



Applying antibiotic disks



Incubation 24h at 37°C




Read the diameter of the inhibition zone



Standardized filter-paper disc-agar diffusion


Procedure



1

→ Transfer at least three to five well-isolated colonies of the same morphological type into nutrient broth tube

2



↓ Incubated between 2 to 6 hrs



حتى اعطيه
هجال
3 للانفصال

4

← Compare the turbidity of the nutrient broth to the 0.5 McFarland standards by either a photometric device or visually.

Standardized filter-paper disc-agar diffusion

Procedure

* بجاوول انا Soap
بمسحكي 100%
surface of area.

5




Dip a sterile cotton swab into a well-mixed saline test



6

streak the entire agar surface horizontally, vertically, and around the outer edge of the plate



Confluent streaking

7

have ← Disk
one antibiotic

Carefully place the provided antibiotic discs onto the plate at equal distances using a sterile forceps and lightly touch each disc to make sure it will stay in place




Standardized filter-paper disc-agar diffusion



Procedure

قاعدة ← شكل Halo zone حول البكتيريا
يعني ان البكتيريا ← Sensitive

• فاستخدمو مسطرة لقياس
Diameter of inhibition
Zone .




Incubation 24h
At 37°C



Standardized filter-paper disc-agar diffusion



Results

Antibiotic	Disk Conc.	Diameter of zone of inhibition (ZOI)		
		Resistant	Intermediate	Susceptible
Amikacin	10 µg	≤11	12-13	≥14
Ampicillin	10 µg	≤11	12-13	≥14
Bacitracin	10 units	≤8	9-11	≥13
Cephalothin	30 µg	≤14	15-17	≥18
Chloramphenicol	30 µg	≤12	13-17	≥18
Clindamycin	2 µg	≤14	15-16	≥17
Erythromycin	15 µg	≤13	14-17	≥18
Gentamicin	10 µg	≤12	13-14	≥15
Kanamycin	30 µg	≤13	14-17	≥18
Lincomycin	2 µg	≤9	10-14	≥15
Methicillin	5 µg	≤9	10-13	≥14
Nalidixic acid	30 µg	≤13	14-18	≥19
Neomycin	30 µg	≤12	13-16	≥17
Nitrofurantoin	0.3 mg	≤14	15-16	≥17
Penicillin				
vs. staphylococci	10 units	≤20	21-28	≥29
vs. other organisms	10 units	≤11	12-21	≥22
Polymyxin	300 units	≤8	9-11	≥12
Streptomycin	10 µg	≤11	12-14	≥15
Sulfonamides	0.3 mg	≤12	13-16	≥17
Tetracycline	30 µg	≤14	15-18	≥19
Vancomycin	30 µg	≤9	10-11	≥12



Standardized filter-paper disc-agar diffusion

Results




* بقيس وبقول مقارنه مع ال standard من الشركة .

Minimum Inhibitory concentration (MIC)
&
Minimum lethal concentration (MLC)

MIC: is the lowest concentration of an antimicrobial that will inhibit the visible growth of a microorganism after overnight incubation

MLC (MBC): Is the lowest concentration of an antibacterial agent required to kill a particular bacterium. It can be determined from broth dilution minimum inhibitory concentration (MIC) tests by subculturing to agar plates that do not contain the test agent.



2

Minimum Inhibitory concentration

1


$\times 10^5$ CFU $\times 10^5$ $\times 10^5$ $\times 10^5$ $\times 10^5$ $\times 10^5$ $\times 10^5$ $\times 10^5$

10^5

10^5

4

• why we use serial dilution?




μg/ml 32 16 8 4 2 1 0.5

Antibiotic serial dilution

To determine how much the least concentration able to prevent bacterial growth by naked eye.

3




MIC

5

Minimum Inhibitory concentration

Minimum Inhibitory concentration



0.5 1 2 4 8 16 µg/ml 32


I can't see bacterial growth by naked eye

Sub-culture to agar medium

ناتج يوم بسوف مين الـ 16
واللي ما اعطيت ←
Growth in the surface of area → MLC



Serial Dilution Susceptibility Testing



Clinical applications for the Qualitative Antimicrobial Susceptibility Testing

MICs can also be used to reduce drug dosage and cost of antimicrobial therapy for very susceptible organisms; therefore, drugs with lower MIC scores are .more effective antimicrobial agents

This is important because populations of bacteria exposed to an insufficient concentration of a particular drug or to a broad-spectrum antibiotic (one designed to inhibit many strains of bacteria) can evolve resistance to these drugs. Therefore, MIC scores aid in improving outcomes for patients and preventing evolution of drug-resistant microbial strains

MIC is used for determining treatment for patients suffering from infections such as sepsis, pneumonia, meningitis, endocarditis or osteomyelitis or managing the treatment of high-risk patients such as those suffering from cystic fibrosis or immunocompromised individuals.

Epsilometer test (E-test)

Used as a substitution for the MIC test

Plastic strips with a predefined gradient of

One antibiotic
One antifungal

One strip per antibiotic

Easy to use


Storage at -20°C

Short shelf life, expensive




Epsilonometer test (E-test)

Elliptical zone of inhibition




MIC 0.25 $\mu\text{g/ml}$

Reading E-tests



Ciprofloxacin



* Antimicrobial Susceptibility Test :

- New antibiotics are continuously being developed
 - Different bacteria acquire new resistant genes to the available antibiotics.
- Determine the antibiotic susceptibility or resistance is required to determine most antibiotic

* Methods of antimicrobial Susceptibility Testing.

Standardized Filter-paper disc agar diffusion (Kirby-Bauer method)


Minimum Inhibitory concentration (MIC)

Epsilometer Test (E-Test)

Qualitative (sensitive or not)

Minimum Lethal Concentration (MLC)
Qualitative + Quantitative

(1) Standardized Filter-paper disc agar diffusion :-




* ① E-coli in Macconkey agar.

② Transfer three to four colony and put them in tube of nutrient broth → Incubator (2-6 hrs)

③ put them in why? حواء طيبا حواء تنفعل

④ 0.5 McFarland → solution. Compare turbidity of the standard with broth after 6 hrs



⑤ Dip swab in bath.

⑥ Streak the entire agar (hinton agar) surface horizontally, around the outer edge

↓ حواء طيبا حواء تنفعل

« Confluent streaking »

⑦ Carefully place the provided antibiotic disk on the plate at equal distances.

* halo zone antibiotic in disk → diffusion in surface of agar if this bacteria is sensitive

↓ the antibiotic will kill it

disk حواء طيبا حواء تنفعل

* Minimum Inhibitory Concentration & Minimum Lethal Concentration

MIC → Lowest concentration antibiotic that will inhibit the visible growth.
- we can determine by naked eye.

MLC → Lowest concentration of an antibacterial agent required to kill a particular bacterium.
- we can determine by subculturing to agar plates that don't contain the test agent.

- ① E-coli in Macconkey agar.
- ② Add 3-4 colony in nutrient growth
- ③ Antibiotic serial dilution in different concentration.
- ④ Add Fixed number of bacteria in a tube.
- ⑤ In the next day → see where the least concentration that I can't see the Growth.

Growth وين اقل تركيز لا نشوف فيه نمو

growth قبل antibiotic *
Concentration

المشاهدة Growth من قبل antibiotic ← tube 3 ←
↳ Sub culture to agar medium

specific No Growth
↓
MLC

③ Epsilon meter test (E-test) :

Used → ~~sub~~ substitution For the MIC test.

uses

↳ plastic strips with a predefined gradient

- One antibiotic ← ^{antib}
- One antiFungal

→ From low concentration to higher.

- One strip per antibiotic
- Storage at 20°C
- Easy to use
- Short shelf life, expensive.

Clinical application

for the

Quantitative

Antimicrobial Susceptibility

Testing

MIC & MFC

● MIC Used to reduce drug dosage and cost of antimicrobial therapy For very susceptible organisms



Drug with lower MIC scores are more effective antimicrob

● MIC scores aid in improving outcomes For patient and preventing evolution of drug-resistant microbial agent.



populations of bacteria exposed to an insufficient concentration of a particular drug (one designed to inhibit many strains)

● MIC is used For determining treatment For patients suffering From infection Like 1. sepsis 2. meningitis 3. pneumonia 4. endocarditis or osteomyelitis or managing the treatment of high-risk patient