


# TYPES OF CULTURE MEDIA

MATHHAR AHMAD ABU MORAD MD

DEPARTMENT OF MICROBIOLOGY AND IMMUNOLOGY

FACULTY OF MEDICINE, MU'TAH UNIVERSITY

## LAB 4



# Purpose

To become familiar with the selective and differential media used to identify the infections associated bacteria


## Principle

- Bacteria and other microbes have particular requirements for growth. Therefore, in order to successfully grow the bacteria in lab so that we can stain and identify them, we must provide an environment that is suitable for growth.
- Growth media are used to cultivate bacteria because it contains essential:
  - ✓ Necessary nutrients
  - ✓ Moisture
  - ✓ pH to support microbial growth

# Streaking Microbial Cultures on Agar Plates

Agar plate streaking are an essential tool in microbiology. They allow bacteria and fungi to grow on a semi-solid surface to produce discrete colonies. These colonies can be used to help identify the organism



## Quadrant Streak



## Plate streaking technique

# Streaking Microbial Cultures on Agar Plates

## Plate streaking technique



## Quadrant Streak

# Overview of bacterial infections

## Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

## Otitis media

- *Streptococcus pneumoniae*

## Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:


- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

## Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*



## Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

## Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

## Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

## Gastritis

- *Helicobacter pylori*

## Food poisoning

- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

## Sexually transmitted diseases

- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

## Urinary tract infections

- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*

## Types of media

```
graph LR; A[Types of media] --- B[Basic media]; A --- C[Enriched media]; A --- D[Enrichment media]; A --- E[Differential media]; A --- F[Selective media]; A --- G[Transport media];
```

**Basic media**

**Enriched media**

**Enrichment media**

**Differential media**

**Selective media**

**Transport media**

# BASIC MEDIA

- Simple media
- Support growth of microorganisms
- No special nutritional requirements
- **Examples:**

Nutrient agar

Nutrient broth




**Nutrient agar**



**Nutrient broth**

# Types of media





# Overview of bacterial infections

## Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

## Otitis media

- *Streptococcus pneumoniae*

## Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:


- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

## Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*



## Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

## Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

## Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

## Gastritis

- *Helicobacter pylori*

## Food poisoning

- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

## Sexually transmitted diseases

- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

## Urinary tract infections

- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*

## Microbiological Analysis of Urine Specimens

<i>Stapylococcus aureus</i>		
<i>Streptococcus pyogenes</i>		
<b>Bacteria</b>	<b>Gram positive</b>	Enterococci
	<b>Gram negative</b>	
		<i>Streptococcus facalis</i> <i>Streptococcus faecium</i>
		<i>Escherichia coli</i> <i>Pseudomonas aeruginosa</i> <i>Proteus vulgaris</i> <i>Klebsiella pneumoniae</i>
<b>Viruses</b> _____	<b>Venereal Disease</b>	<i>Treponema pallidum</i> <i>Neisseria gonorrhoeae</i> <i>Hemophilus ducreyi</i>
	—	<i>Calymnatobacterium granulomatis</i> <i>Herpes hominus (type 11)</i>
<b>Fungi</b>	<i>Candida albicans</i> <i>Blastomyces dermatitidis</i> <i>Coccidioides bancrofti</i>	
<b>Protozoa</b>	<i>Trichomonas vaginalis</i> <i>Entameoba histolytica</i>	

# Urine analysis

Midstream urine sample  
Sample inoculation

Blood agar

Significant growth

Gram stain

Gram positive cocci

catalase

positive

Staphylococci

Mannitol salt agar

negative

Streptococci

MacConkey agar

Abundant growth


Gram stain

Gram negative bacilli

*Escherichia coli, Pseudomonas aeruginosa*  
*Proteus vulgaris, Klebsiella pneumoniae*

Biochemical reactions

# MacConkey agar is a selective and differential media for Enterobacteriaceae






↓

*Escherichia coli*  
*Klebsiella spp*  
*Enterobacter spp*  
*Citrobacter spp*

↓

*Salmonella spp*  
*Shigella spp*  
*Proteus spp*  
*Yersinia spp*

# MacConkey agar



# Blood agar (BA)

Enriched medium: containing peptones, yeast extracts, liver or heart extracts (depending on the medium), and blood.



Some bacteria produce an enzyme called hemolysin that is able to lyse RBCs (hemolysis)



Differential medium: containing blood

If hemolysin is produced by the bacteria it will be secreted into the medium and the RBCs will be lysed

Growth on BA differentiates between the three groups of Bacteria:

- 1- Alpha hemolytic bacteria
- 2- Gamma hemolytic bacteria
- 3- Beta hemolytic bacteria

# Beta hemolysis = Complete hemolysis



# Alpha hemolysis

Hemoglobin containing  
 $\text{Fe}^{2+}$  (ferrous)




hydrogen peroxide  
produced by the bacterium



**Oxidation of  $\text{Fe}^{2+}$   
into  $\text{Fe}^{3+}$  (ferric) state**



Hemoglobin converted into  
methemoglobin (greenish color)




Incomplete (partial) lysis of RBCs



# Gamma hemolysis

No hemolysis, and no change in the medium




The three types of hemolysis



# The hemolytic pattern of different Streptococci

## Streptococcus



# Mannitol salt agar

Selective and Differential for Staphylococci





- **Selective agent:** 7.5% NaCl
- **Differential agent:** mannitol to differentiate between mannitol Fermenters and non-fermenters
- pH indicator: Phenol red

# Mannitol salt agar

MSA



Non-cultured



Blood agar



Beta hemolytic Staphylococci

# Overview of bacterial infections

## Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

## Otitis media

- *Streptococcus pneumoniae*

## Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:


- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

## Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*



## Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

## Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

## Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

## Gastritis

- *Helicobacter pylori*

## Food poisoning

- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

## Sexually transmitted diseases

- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

## Urinary tract infections


- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*

# Processing of stool samples

Patient with diarrhea



Stool sample received in the lab



Culture for bacteria



Microscopy for parasites

# *Salmonella -Shigella* agar (SS agar)

## **Purpose**

For isolation and differentiation of *Salmonella* & *Shigella*

## **Components**

- ✓ the presence of **brilliant green, sodium citrate and bile salts** which completely inhibit the growth of Gram-positive bacteria and partially inhibit the growth of Enterobacteriaceae and Proteus.
- ✓ **Lactose**: carbon source
- ✓ **Neutral red**: pH indicator, red in acidic conditions

# *Salmonella Shigella* agar (SS agar)

Why black colonies?

Due to the production of FeS (ferrous sulfide forming black precipitate presented by black-centered colonies)

SS agar




- **Sodium thiosulfate** ( $\text{Na}_2\text{S}_2\text{O}_3$ ): sulfur source
- **Fe<sup>3+</sup>** (ferric) H<sub>2</sub>S indicator



$\text{Na}_2\text{S}_2\text{O}_3$  + thiosulfate reductase  $\longrightarrow$  sulfite + H<sub>2</sub>S

H<sub>2</sub>S + Fe<sup>3+</sup>  $\longrightarrow$  FeS (black precipitate presented by black-centered colonies)





# *Salmonella Shigella* agar (SS agar)

## Results

- ✓ Lactose fermenters: pink to red colonies (few can grow)
- ✓ Non lactose fermenters: translucent, colorless colonies with or without black centers



***Shigella***: colorless colonies without black centers

**Lactose fermenter flora:**  
pink to red colonies



***Salmonella***:  
colorless colonies  
with black centers



# Overview of bacterial infections

## Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

## Otitis media

- *Streptococcus pneumoniae*

## Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:


- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

## Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*



## Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

## Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

## Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

## Gastritis

- *Helicobacter pylori*

## Food poisoning

- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

## Sexually transmitted diseases

- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

## Urinary tract infections

- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*


# Cholera identification

## Identification

- ✓ Thiosulfate citrate bile salt sucrose agar or TCBS agar
- ✓ The medium is alkaline which enhances the growth of *Vibrio* species

## Important components

- ✓ Sucrose: sugar source
- ✓ Bromothymol blue: pH indicator
  - pH < 6.0 - yellow
  - pH > 7.6 - blue



# *Cholera* identification

## Results


- ✓ *Vibrio cholera*: Ferment sucrose smooth yellow colonies
- ✓ *Vibrio parahemolyticus*: non-sucrose fermenter, green colonies



TCBS media




*V. cholera*



*V. parahemolyticus*

# *Cholera* identification



# Overview of bacterial infections

## Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

## Otitis media

- *Streptococcus pneumoniae*

## Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:


- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

## Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*



## Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

## Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

## Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

## Gastritis

- *Helicobacter pylori*

## Food poisoning

- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

## Sexually transmitted diseases



- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

## Urinary tract infections

- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*

# Löwenstein–Jensen (LJ) medium

- Is a growth medium specially used for culture of *Mycobacterium*, notably *Mycobacterium tuberculosis*.



M.tuberculosis produces rough and tough colonies

# Löwenstein–Jensen (LJ) medium


**Penicillin** and **Nalidixic acid** (FD053) along with **malachite green** prevents growth of the majority of contaminants surviving decontamination of the specimen while encouraging earliest possible growth of Mycobacteria

Composition	Ingredients Gms / 600 ml
L-Asparagine	<b>3.600</b>
Monopotassium phosphate	<b>2.400</b>
Magnesium sulphate	<b>0.240</b>
Magnesium citrate	<b>0.600</b>
Potato starch, soluble	<b>30.000</b>
Malachite green	<b>0.40</b>




## Cultivation of Anaerobic Bacteria – Thioglycollate medium

**Thioglycollate** is a multipurpose, enriched, differential medium used primarily to determine the oxygen requirements of microorganisms. It acts as an enrichment broth which is most frequently used in diagnostic bacteriology. This broth supports the growth of anaerobes, aerobes, facultative anaerobes, microaerophilic, and aerotolerant microorganisms.



# Anaerobic jars




Anaerobic candle jar

# Löffler's medium

Is a special substance used to grow *Corynebacterium diphtheriae* bacilli to confirm the diagnosis.



Gram-positive rod-shaped bacteria that are straight or slightly curved. The bacteria group together in a characteristic way (Chinese letters)



# Overview of bacterial infections

## Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

## Otitis media

- *Streptococcus pneumoniae*

## Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:


- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

## Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*



## Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

## Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

## Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

## Gastritis

- *Helicobacter pylori*

## Food poisoning

- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*


## Sexually transmitted diseases

- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

## Urinary tract infections

- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*

# CSF analysis



# Chocolate agar

- Used to isolate *Haemophilus influenzae*
- Is a hemolysed blood either by heating blood to 80°C or using enzyme treatment
- Treatment result in browning of the medium, therefore, it is called chocolate agar.



Chocolate agar



Haemophilus influenzae growth on Chocolate agar



**Thank you**