


Microbiology

Lab 4



- Blood agar \Rightarrow Gram + cocci
- Manitol salt sugar \Rightarrow staph
- MacConkey agar \Rightarrow Gram - bacilli
- SS agar \Rightarrow Salmonella - Shigella.
- TCBS agar \Rightarrow Cholera .
- LJ medium \Rightarrow mycobacterium
- Chocolate agar \Rightarrow Haemophilus influenzae.
- Anaerobic jars \Rightarrow anaerobic bacteria



General Microbiology Lab

Types of Culture Media 2021-2022

Lab 4

Dr. Mohammad Odaibat
Department of Microbiology and Pathology
Faculty of Medicine, Mutah University

Purpose

To become familiar with the selective and differential media used to identify the infections associated bacteria

Principle

• Bacteria and other microbes have particular requirements for growth. Therefore, in order to successfully grow the bacteria in lab so that we can stain and identify them, we must provide an environment that is suitable for growth.

- Growth media are used to cultivate bacteria because it contains essential:
 - ✓ Necessary nutrients
 - ✓ Moisture
 - ✓ pH to support microbial growth

Not all bacteria need the same growth requirement

Streaking Microbial Cultures on Agar Plates

Agar plate streaking are an essential tool in microbiology. They allow bacteria and fungi to grow on a semi-solid surface to produce discrete colonies. These colonies can be used to help identify the organism

Quadrant Streak





Plate streaking technique

To know whether the infection is caused by one type of bacteria or more


Mono infection (one type) /// Mixed infection (two colonies)

Streaking Microbial Cultures On Agar Plates

Plate streaking technique



Quadrant Streak



Overview of bacterial infections

Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

Otitis media

- *Streptococcus pneumoniae*

Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:

- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*

Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

Gastritis

- *Helicobacter pylori*

Food poisoning


- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

Sexually transmitted diseases

- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

Urinary tract infections

- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*



Types of media

Types of media

Enriched media

مزودة بفيتامينات ومواد غذائية كثيرة

معظم أنواع البكتيريا

contains **specific growth factors** needed by **fastidious bacteria** to support their growth.

Examples

-blood agar

-chocolate agar

Selective media

تختار نوع واحد من البكتيريا

used to **select (isolate)** **specific group of bacteria**

Differential media

تميز
Species

these can distinguish among morphologically and biochemically related groups of organisms.

Overview of bacterial infections

Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

Otitis media

- *Streptococcus pneumoniae*

Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:


- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*



Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

Gastritis

- *Helicobacter pylori*

Food poisoning

- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

Sexually transmitted diseases

- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

Urinary tract infections

- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*

UTI

Most common cause of UTI /80%

Microbiological Analysis of Urine Specimens

| | | |
|-----------------|---|--|
| Bacteria | Gram positive | <i>Stapylococcus aureus</i> <i>Streptococcus pyogenes</i> <i>Streptococcus facalis</i> <i>Streptococcus faecium</i> |
| | Gram negative | <i>Escherichia coli</i> <i>Pseudomonas aeruginosa</i> <i>Proteus vulgaris</i> <i>Klebsiella pneumoniae</i> |
| Viruses | Venereal Disease | <i>Treponema pallidum</i> <i>Neisseria gonorrhoeae</i> <i>Hemophilus ducreyi</i> |
| | — | <i>Herpes hominus (type 11)</i> |
| Fungi | <i>Candida albicans</i> <i>Blastomyces dermatitidis</i> <i>Coccidioides bancrofti</i> | |
| Protozoa | <i>Trichomonas vaginalis</i> <i>Entameoba histolytica</i> | |

Urine analysis

Midstream urine sample
Sample inoculation

The second third of the urine. No contamination with the bacteria in urethrae

Blood agar

MacConkey agar

Selective and differential media

Non selective/enriched media

Significant growth

Abundant growth

Gram stain

Gram stain

ثاني يوم

Gram positive cocci

Gram negative bacilli

catalase

Escherichia coli, Pseudomonas aeruginosa
Proteus vulgaris, Klebsiella pneumoniae

positive

negative

Biochemical reactions

بعد اسبوع

Staphylococci


Streptococci

بتحدد المسبب لل infection

Mannitol salt agar

Differential Media for staph

MacConkey agar is a selective and differential media for Enterobacteriaceae



Gram negative lactose fermenter



Not enough!
We need

Biochemical reaction

- Escherichia coli*
- Klebsiella spp*
- Enterobacter spp*
- Citrobacter spp*

- Salmonella spp*
- Shigella spp*
- Proteus spp*
- Yersinia spp*

MacConkey agar



Blood agar (BA)

Enriched medium: containing peptones, yeast extracts, liver or heart extracts (depending on the medium), and blood.



Some bacteria produce an enzyme called hemolysin that is able to lyse RBCs (hemolysis)

Differential medium: containing blood

If hemolysin is produced by the bacteria it will be secreted into the medium and the RBCs will be lysed

Growth on BA differentiates between the three groups of Bacteria:

- 1- Alpha hemolytic bacteria Partial hemolysis
- 2- Gamma hemolytic bacteria Complete hemolysis
- 3- Beta hemolytic bacteria No hemolysis


Notes in the previous slide

*A growth medium is considered **differential** if, when specific microbes are present, the medium or bacterial colonies themselves exhibit a color change that provides **information** about their identity.

*Blood agar (**BAP**) is a **differential** growth medium which microbiologists use to distinguish clinically significant bacteria from **throat** and **sputum** cultures. BAP contains **5% sheep blood**. Certain bacteria produce exotoxins called hemolysins, which act on the red blood cells to lyse, or break them down.

Microbiologist use differential media to identify and isolate specific bacteria. An example of this is the bacteria *Streptococcus pyogenes*, which causes strep throat. You can grow these bacteria on a complex media such as nutrient agar, but if other bacteria are also growing on that agar, it is very difficult to distinguish one bacterial colony from another without the use of microscopic examination and special staining techniques. If you grow it on blood agar, though, it will destroy the red blood cells in a process called beta-hemolysis, and other cells will not cause this reaction, which makes identifying *Streptococcus pyogenes* much easier.

Beta hemolysis = Complete hemolysis



* Transparent hallos around bacterial growth

Alpha hemolysis



Greenish

Hemoglobin containing
 Fe^{2+} (ferrous)




hydrogen peroxide
produced by the bacterium



**Oxidation of Fe^{2+}
into Fe^{3+} (ferric) state**



Hemoglobin converted into
methemoglobin (greenish color)




Incomplete (partial) lysis of RBCs

Methemoglobin (pronounced "met-hemoglobin") is a form of the oxygen-carrying metalloprotein hemoglobin, in which the iron in the heme group is in the Fe^{3+} (ferric) state, not the Fe^{2+} (ferrous) of normal hemoglobin. Methemoglobin cannot bind oxygen, unlike oxyhemoglobin. It is bluish chocolate-brown in color. In human blood a trace amount of methemoglobin is normally produced spontaneously. But when it is present in excess the blood becomes abnormally dark bluish brown. The NADH-dependent enzyme methemoglobin reductase (diaphorase I) is responsible for converting methemoglobin back to hemoglobin.

Gamma hemolysis


No hemolysis, and no change in the medium



The three types of hemolysis

The hemolytic pattern of different Streptococci

Streptococcus



Urine analysis

Midstream urine sample

Sample inoculation

Blood agar

Significant growth

Gram stain

Gram positive cocci

catalase

positive

Staphylococci

Mannitol salt agar

MacConkey agar

Abundant growth

Gram stain

Gram negative bacilli

Escherichia coli, Pseudomonas aeruginosa
Proteus vulgaris, Klebsiella pneumoniae

Biochemical reactions

Mannitol salt agar

Selective and Differential for Staphylococci



- **Selective agent: 7.5% NaCl**
- **Differential agent:** mannitol to differentiate between mannitol Fermenters and nonfermenters
- **pH indicator: Phenol red**

Staphylococcus aureus

Other staphylococcus


Mannitol salt agar

Halophilic

MSA



Non-cultured



Blood agar



Beta hemolytic Staphylococci

Overview of bacterial infections

Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

Otitis media

- *Streptococcus pneumoniae*

Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:

- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*

Sexually transmitted diseases

- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

Gastritis


- *Helicobacter pylori*

Food poisoning

- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

Urinary tract infections



- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*



Processing of stool samples

Patient with diarrhea

Stool sample received in the lab



Culture for bacteria



Microscopy for parasites

قسط
Parasitology

Salmonella -Shigella agar (SS agar)

Purpose

For isolation and differentiation of *Salmonella* & *Shigella*

Components

- ✓ Bile salt: inhibit the growth of gram positive bacteria (selective agent)
- ✓ Lactose: carbon source
- ✓ Neutral red: pH indicator, red in acidic conditions


Salmonella Shigella agar (SS agar)

Why black colonies?
Salmonella

Due to the production of **FeS** (ferrous sulfide forming black precipitate presented by black-centered colonies)

SS agar

- Sodium thiosulfate ($\text{Na}_2\text{S}_2\text{O}_3$): sulfur source
- Fe^{3+} (ferric) H_2S indicator



نستغل وجود هذا الانزيم في ال salmonella حتى نكتشفها

Salmonella Shigella agar (SS agar)

Results

- ✓ **Lactose fermenters:** pink to red colonies (few can grow)
- ✓ **Non lactose fermenters:** translucent, colorless colonies with or without black centers



Shigella: colorless colonies without black centers

Lactose fermenter flora:
pink to red colonies



Salmonella:
colorless colonies with black centers

Owl eye



Overview of bacterial infections

Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

Otitis media

- *Streptococcus pneumoniae*

Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:


- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*



Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

Gastritis

- *Helicobacter pylori*

Food poisoning

- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

Sexually transmitted diseases

- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

Urinary tract infections

- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*

Cholera identification

Identification

- ✓ Thiosulfate citrate bile salt sucrose agar or TCBS agar
- ✓ The medium is alkaline (pH 8.6) which enhances the growth of *Vibrio* species

Important components


- ✓ Sucrose: sugar source
- ✓ Bromothymol blue: pH indicator
 - pH < 6.0 - yellow
 - pH > 7.6 - blue

Cholera identification

Results

- ✓ *Vibrio cholera*: Ferment sucrose smooth yellow colonies
- ✓ *Vibrio parahemolyticus*: non-sucrose fermenter, green colonies

لون ال media



TCBS media



V. cholera



V. parahemolyticus

Overview of bacterial infections

Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

Otitis media

- *Streptococcus pneumoniae*

Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:


- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*



Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

Gastritis

- *Helicobacter pylori*

Food poisoning

- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

Sexually transmitted diseases


- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

Urinary tract infections

- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*

Löwenstein–Jensen (LJ) medium

- Is a growth medium specially used for culture of *Mycobacterium*, notably *Mycobacterium tuberculosis*.



M.tuberculosis produces rough and tough colonies
Because of the micolic acid

Penicillin and Nalidixic acid along with malachite green prevents growth of the majority of contaminants surviving decontamination of the specimen while encouraging earliest possible growth of Mycobacteria

Overview of bacterial infections

Bacterial meningitis

- *Streptococcus pneumoniae*
- *Neisseria meningitidis*
- *Haemophilus influenzae*
- *Streptococcus agalactiae*
- *Listeria monocytogenes*

Chocolate agar

يعيشوا على مواد داخل ال RBCs

Otitis media

- *Streptococcus pneumoniae*

Pneumonia

Community-acquired:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*

Atypical:

- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*
- *Legionella pneumophila*

Tuberculosis

- *Mycobacterium tuberculosis*

Skin infections

- *Staphylococcus aureus*
- *Streptococcus pyogenes*
- *Pseudomonas aeruginosa*

Sexually transmitted diseases

- *Chlamydia trachomatis*
- *Neisseria gonorrhoeae*
- *Treponema pallidum*
- *Ureaplasma urealyticum*
- *Haemophilus ducreyi*

Eye infections

- *Staphylococcus aureus*
- *Neisseria gonorrhoeae*
- *Chlamydia trachomatis*

Sinusitis

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*

Upper respiratory tract infection

- *Streptococcus pyogenes*
- *Haemophilus influenzae*

Gastritis


- *Helicobacter pylori*

Food poisoning


- *Campylobacter jejuni*
- *Salmonella*
- *Shigella*
- *Clostridium*
- *Staphylococcus aureus*
- *Escherichia coli*

Urinary tract infections

- *Escherichia coli*
- Other Enterobacteriaceae
- *Staphylococcus saprophyticus*
- *Pseudomonas aeruginosa*



CSF analysis



Chocolate agar

- Used to isolate *Haemophilus influenzae*
- Is a hemolysed blood either by heating blood to 80°C or using enzyme treatment
- Treatment result in browning of the medium, therefore, it is called chocolate agar.



Chocolate agar




Haemophilus influenzae growth on Chocolate agar

Hemophilus requires both X factor (hematin, containing iron) and V factor (NAD). These should be released from blood cells and therefore we use heating. Heating is also important to inactivate NADase

Anaerobic gars

بنزرع فيه ال
Anaerobic bacteria



Anaerobic candle jar